


Product Datasheet

RNA Polymerase II, Phospho (S2) Antibody (orb1523052)


Description

RNA Polymerase II, Phospho (S2) Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Surrounding serine 2
Target	RNA Polymerase II, Phospho (S2)
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P24928
NCBI	NP_000928.1
Dilution Range	WB - 1:500 - 1:2,500; IP - 2 - 10 µg/mg lysate; IHC - 1:2,000 - 1:10,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human Phospho RNA Polymeras...



Detection of human Phospho RNA Polymeras...



Detection of mouse Phospho RNA Polymeras...