


Product Datasheet

HRP-2 Antibody (orb1523225)

Description	HRP-2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Between 621 and C-terminus
Target	HRP-2
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q7Z4V5
NCBI	NP_001001520.1
Dilution Range	WB - 1:500 - 1:2,500; IP - Not recommended
Expiration Date	12 months from date of receipt.



Detection of human HRP-2 by western blot...