



Product Datasheet

TRMT1 Antibody (orb1523429)

Description	TRMT1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 609 and C-terminus
Target	TRMT1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9NXH9
NCBI	NP_060192.1
Dilution Range	WB - 1:500 - 1:2,500; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human TRMT1 by western blot...



Detection of human TRMT1 by western blot...