


Product Datasheet

MMS21 Antibody (orb1523559)

Description	MMS21 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 197 and C-terminus
Target	MMS21
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q96MF7
NCBI	NP_775956.1
Dilution Range	WB - Not recommended. A304-129A has not performed satisfactorily when used for WB of MMS21 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of MMS21.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human MMS21 by western blot...