


Product Datasheet

RNF216 Antibody (orb1523591)

Description	RNF216 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 775 and 825
Target	RNF216
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9NWF9
NCBI	NP_996999.1
Dilution Range	WB - Not recommended. Use rabbit anti-RNF216 antibody A304-111A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human RNF216 by western blo...