



Product Datasheet

MKI67IP Antibody (orb1523699)

Description	MKI67IP Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 243 and C-terminus
Target	MKI67IP
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9BYG3
NCBI	NP_115766.3
Dilution Range	WB - 1:5,000 - 1:25,000; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human MKI67IP by western bl...



Detection of human MKI67IP by western bl...