


Product Datasheet

HO-1 Antibody (orb1524424)

Description	HO-1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 1 and 50
Target	HO-1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P09601
NCBI	NP_002124.1
Dilution Range	WB - Not recommended; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human HO-1 by western blot ...