


Product Datasheet

ERK1/2, Phospho (T202/Y204, T185/Y187) Antibody (orb1524515)

Description	ERK1/2, Phospho (T202/Y204, T185/Y187) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Surrounding threonine 185 and tyrosine 187
Target	ERK1/2, Phospho (T202/Y204, T185/Y187)
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P27361 , P28482
NCBI	NP_002737.2
Dilution Range	WB - 1:1,000 - 1:5,000; IP - Not recommended
Expiration Date	12 months from date of receipt.



Detection of human phospho-ERK1 and ERK2...