


Product Datasheet

FYCO1 Antibody (orb1525982)

Description	FYCO1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 325 and 375
Target	FYCO1
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9BQS8
NCBI	NP_078789.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 5 - 15 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human FYCO1 by western blot...