


Product Datasheet

TNIK Antibody (orb1526159)

Description	TNIK Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 775 and 825
Target	TNIK
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9UKE5
NCBI	NP_055843.1
Dilution Range	WB - Not recommended. Use rabbit anti-TNIK antibody A302-695A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human TNIK by western blot ...