


Product Datasheet

MERIT40 Antibody (orb1526484)

Description	MERIT40 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 1 and 50
Target	MERIT40
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9NWX8
NCBI	NP_054892.2
Dilution Range	WB - Not recommended. orb1526484 has not performed satisfactorily when used for WB of MERIT40 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of MERIT40.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human MERIT40 by western bl...