



Product Datasheet

BRD1 Antibody (orb1526754)

| | |
|----------------------------|--|
| Description | BRD1 Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IP, WB |
| Immunogen | Between 775 and 825 |
| Target | BRD1 |
| Preservatives | Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide |
| Form/Appearance | Whole IgG |
| Concentration | 200 µg/ml |
| Storage | 2 - 8°C |
| Note | For research use only |
| Application notes | Format: Whole IgG |
| Clonality | Polyclonal |
| Uniprot ID | O95696 |
| NCBI | NP_055392.1 |
| Dilution Range | WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate |
| Expiration Date | 12 months from date of receipt. |



Detection of human BRD1 by western blot....



Detection of human BRD1 by western blot ...