


Product Datasheet

PPP1R14B Antibody (orb1526755)

Description	PPP1R14B Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 97 and C-terminus
Target	PPP1R14B
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q96C90
NCBI	NP_619634.1
Dilution Range	WB - Not recommended. A302-365A has not performed satisfactorily when used for WB of PPP1R14B in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of PPP1R14B.; IP - 5 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human PPP1R14B by western b...