



Product Datasheet

TIGAR Antibody (orb1526893)

| | |
|----------------------------|---|
| Description | TIGAR Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC, IP |
| Immunogen | Between 220 and C-terminus |
| Target | TIGAR |
| Preservatives | Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide |
| Form/Appearance | Whole IgG |
| Concentration | 1000 µg/ml |
| Storage | 2 - 8°C |
| Note | For research use only |
| Application notes | Format: Whole IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q9NQ88 |
| NCBI | NP_065108.1 |
| Dilution Range | WB - Not recommended. Use rabbit anti-TIGAR antibody A302-285A.; IP - 10 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections. |
| Expiration Date | 12 months from date of receipt. |



Detection of human TIGAR by western blot...



Detection of human TIGAR by immunohistoc...