


Product Datasheet

Abraxas Antibody (orb1527072)

Description	Abraxas Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 359 and C-terminus
Target	Abraxas
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q6UWZ7
NCBI	NP_620775.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 5 - 15 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human Abraxas by western blotting.