



## Product Datasheet

# ZAP70, Phospho (Y493/Y494) Antibody (orb1527393)

<b>Description</b>	ZAP70, Phospho (Y493/Y494) Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IP
<b>Immunogen</b>	Surrounding tyrosine 493/494
<b>Target</b>	ZAP70, Phospho (Y493/Y494)
<b>Preservatives</b>	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
<b>Form/Appearance</b>	Whole IgG
<b>Concentration</b>	1000 µg/ml
<b>Storage</b>	2 - 8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Format: Whole IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P43403</a>
<b>NCBI</b>	<a href="#">NP_001070.2</a>
<b>Dilution Range</b>	WB - Not recommended. orb1527393 has not performed satisfactorily when used for WB of phospho ZAP70 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of ZAP70.; IP - 2 - 5 µg/mg lysate
<b>Expiration Date</b>	12 months from date of receipt.



Detection of Phosphorylation of human ZA...



Detection of Phosphorylation of human ZA...