



Product Datasheet

ZAP70, Phospho (Y319) Antibody (orb1527397)

Description	ZAP70, Phospho (Y319) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	Surrounding tyrosine 319
Target	ZAP70, Phospho (Y319)
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P43403
NCBI	NP_001070.2
Dilution Range	WB - Not recommended. orb1527397 has not performed satisfactorily when used for WB of phospho ZAP70 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of ZAP70.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of Phosphorylation of human ZA...



Detection of Phosphorylation of human ZA...