



## Product Datasheet

### ZAK Antibody (orb1527418)

<b>Description</b>	ZAK Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, IP
<b>Immunogen</b>	Between 750 and C-terminus
<b>Target</b>	ZAK
<b>Preservatives</b>	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
<b>Form/Appearance</b>	Whole IgG
<b>Concentration</b>	1000 µg/ml
<b>Storage</b>	2 - 8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Format: Whole IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q9NYL2</a>
<b>NCBI</b>	<a href="#">NP_057737.2</a>
<b>Dilution Range</b>	WB - Not recommended. Use rabbit anti-ZAK antibody A301-993A.; IP - 2 - 5 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
<b>Expiration Date</b>	12 months from date of receipt.



Detection of human ZAK by western blot o...



Detection of human ZAK by immunohistoche...