



Product Datasheet

FANCA Antibody (orb1527443)


| | |
|----------------------------|--|
| Description | FANCA Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | FC, IP, WB |
| Immunogen | Between 1405 and C-terminus |
| Target | FANCA |
| Preservatives | Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide |
| Form/Appearance | Whole IgG |
| Concentration | 200 µg/ml |
| Storage | 2 - 8°C |
| Note | For research use only |
| Application notes | Format: Whole IgG |
| Clonality | Polyclonal |
| Uniprot ID | O15360 |
| NCBI | NP_000126.2 |
| Dilution Range | WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate |
| Expiration Date | 12 months from date of receipt. |



Detection of human FANCA by western blot...



Detection of human FANCA by western blot...



Detection of human FANCA (shaded) in HDL...