



Product Datasheet

GRLF1/p190RhoGAP Antibody (orb1527883)

Description	GRLF1/p190RhoGAP Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 1200 and 1250
Target	GRLF1/p190RhoGAP
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9NRY4
NCBI	EAW57452.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of mouse GRLF1/p190RhoGAP by w...



Detection of human GRLF1/p190RhoGAP by w...