


Product Datasheet

NUMBL Antibody (orb1527914)

| | |
|----------------------------|--|
| Description | NUMBL Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IP |
| Immunogen | Between 559 and C-terminus |
| Target | NUMBL |
| Preservatives | Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide |
| Form/Appearance | Whole IgG |
| Concentration | 1000 µg/ml |
| Storage | 2 - 8°C |
| Note | For research use only |
| Application notes | Format: Whole IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q9Y6R0 |
| NCBI | NP_004747.1 |
| Dilution Range | WB - Not recommended. A301-719A has not performed satisfactorily when used for WB of NUMBL in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of NUMBL.; IP - 2 - 5 µg/mg lysate |
| Expiration Date | 12 months from date of receipt. |



Detection of human NUMBL1 by western blo...