


Product Datasheet


ENH/PDLIM5 Antibody (orb1527939)

Description	ENH/PDLIM5 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 546 and C-terminus
Target	ENH/PDLIM5
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q96HC4
NCBI	NP_006448.2
Dilution Range	WB - 1:1,000 - 1:2,000; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human and mouse ENH/PDLIM5

...



Detection of human ENH/PDLIM5 by western...