


Product Datasheet

ENH/PDLIM5 Antibody (orb1527941)

Description	ENH/PDLIM5 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 375 and 425
Target	ENH/PDLIM5
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q96HC4
NCBI	NP_006448.2
Dilution Range	WB - Not recommended. Use rabbit anti-ENH/PDLIM5 antibody A301-704A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human ENH/PDLIM5 by western...