



Product Datasheet

PIP1/WDR84 Antibody (orb1528199)

Description	PIP1/WDR84 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 325 and 375
Target	PIP1/WDR84
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9NWT1
NCBI	NP_060376.2
Dilution Range	WB - 1:1,000 - 1:5,000; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human PIP1/WDR84 by western...



Detection of human PIP1/WDR84 by western...