



Product Datasheet

TRIM33/TIF1gamma Antibody (orb1529123)

Description	TRIM33/TIF1gamma Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Between 1 and 50
Target	TRIM33/TIF1gamma
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9UPN9
NCBI	NP_056990.3
Dilution Range	WB - 1:2,000 - 1:10,000; IP - Not recommended. Use rabbit anti-TRIM33/TIF1gamma antibody A301-059A or A301-060A.
Expiration Date	12 months from date of receipt.



Detection of human TRIM33/TIF1gamma by w...



Detection of human TRIM33/TIF1gamma by w...