



Product Datasheet

Pescadillo Antibody (orb1529421)

Description	Pescadillo Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	Between 25 and 75
Target	Pescadillo
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O00541
NCBI	NP_055118.1
Dilution Range	WB - Not recommended. Use rabbit anti-Pescadillo antibody A300-903A.; IP - 2 - 5 µg/mg lysate; IHC - 1:1,000 - 1:5,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human Pescadillo by immunohistochemistry (IHC).



Detection of human Pescadillo by western blot analysis (WB).