


Product Datasheet

HURP/DLG7 Antibody (orb1529515)

Description	HURP/DLG7 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 150 and 200
Target	HURP/DLG7
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q15398
NCBI	NP_055565.2
Dilution Range	WB - Not recommended. Use rabbit anti-HURP/DLG7 antibody A300-852A or A300-853A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human HURP/DLG7 by western ...