



Product Datasheet

TIF1 Alpha/TRIM24 Antibody (orb1529581)

Description	TIF1 Alpha/TRIM24 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 1000 and C-terminus
Target	TIF1 Alpha/TRIM24
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O15164
NCBI	NP_056989.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human TIF1 alpha/TRIM24 by ...



Detection of human TIF1 alpha/TRIM24 by ...