


Product Datasheet

MAVS/VISA Antibody (orb1529643)

Description	MAVS/VISA Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 275 and 325
Target	MAVS/VISA
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q7Z434
NCBI	AAH44952.1
Dilution Range	WB - Not recommended. Use rabbit anti-MAVS antibody A300-782A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human MAVS/VISA by western ...