


Product Datasheet

RPS6, Phospho (S235/S236) Antibody (orb1529954)

Description	RPS6, Phospho (S235/S236) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Surrounding serine 235/236
Target	RPS6, Phospho (S235/S236)
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P62753
NCBI	NP_001001.2
Dilution Range	WB - 1:2,500 - 1:10,000
Expiration Date	12 months from date of receipt.



Detection of human Phospho RPS6 (S235/S236)