


Product Datasheet

MST3/STK24 Antibody (orb1529987)

Description	MST3/STK24 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 325 and 375
Target	MST3/STK24
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9Y6E0
NCBI	NP_003567.2
Dilution Range	WB - 1:1,000 - 1:10,000; IP - 1 - 4 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human MST3(STK24) by WB and...