



Product Datasheet

BRIP1 Antibody (orb1529989)

Description	BRIP1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 1199 and C-terminus
Target	BRIP1
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9BX63
NCBI	NP_114432.1
Dilution Range	WB - 1:1,000 - 1:10,000; IP - 1 - 4 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human BRIP1/BACH1 by wester...



Detection of human BRIP1/BACH1 by wester...