


Product Datasheet

SMG1 Antibody (orb1530265)

Description	SMG1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Between 1625 and 1675
Target	SMG1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q96Q15
NCBI	Q96QV0
Dilution Range	WB - 1:2,000 - 1:10,000; IP - Not recommended. A300-394A is recommended.
Expiration Date	12 months from date of receipt.



Detection of human SMG1 by western blot....