



Product Datasheet

DHX15 Antibody (orb1530268)

Description	DHX15 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 745 and C-terminus
Target	DHX15
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O43143
NCBI	NP_001349.2
Dilution Range	WB - 1:5,000 - 1:15,000; IP - 1 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human DHX15 by western blot...



Detection of human DHX15 by WB and IP.Sa...