



Product Datasheet

TSC1 Antibody (orb1530367)


Description	TSC1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 1114 and C-terminus
Target	TSC1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q92574
NCBI	NP_000359.1
Dilution Range	WB - 1:5,000 - 1:25,000; IP - 1 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human and mouse TSC1 by wes...



Detection of human and mouse TSC1 by wes...



Detection of human TSC1 by western blot ...