



Product Datasheet

KIAA0082 Antibody (orb1530381)

Description	KIAA0082 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 785 and C-terminus
Target	KIAA0082
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8N1G2
NCBI	NP_055865.1
Dilution Range	WB - 1:5,000 - 1:25,000; IP - 1 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human KIAA0082 by western b...



Detection of human KIAA0082 by WB and IP...