



Product Datasheet

PCNA Antibody (orb1530429)


Description	PCNA Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Between 211 and C-terminus
Target	PCNA
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P12004
NCBI	NP_872590.1
Dilution Range	WB - 1:5,000 - 1:20,000; IP - 2 - 10 µg/mg lysate; IHC - 1:2,000 to 1:10,000. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.; IHC-IF - 1:1,000 - 1:5,000
Expiration Date	12 months from date of receipt.



Detection of human PCNA by western blot....



Detection of mouse PCNA by western blot....



Detection of mouse PCNA by immunohistochem...