


Product Datasheet

WAPL Antibody (orb1530446)

Description	WAPL Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Between 1140 and C-terminus
Target	WAPL
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q7Z5K2
NCBI	NP_055860.1
Dilution Range	WB - 1:10,000 - 1:25,000
Expiration Date	12 months from date of receipt.



Detection of human WAPL by western blot....