


Product Datasheet

SAH3 Antibody (orb1530457)

Description	SAH3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 1 and 50
Target	SAH3
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q96HN2
NCBI	NP_056143.1
Dilution Range	WB - 1:1,000 - 1:10,000; IP - 1 - 4 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human SAH3 by WB and IP.Sam...