


Product Datasheet

WASP, Phospho (S483/S484) Antibody (orb1530544)

Description	WASP, Phospho (S483/S484) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Surrounding serine 483/484
Target	WASP, Phospho (S483/S484)
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P42768
NCBI	NP_000368.1
Dilution Range	WB - 1:1,000 - 1:10,000
Expiration Date	12 months from date of receipt.



Detection of human Phospho WASP (S483/S484)