



Product Datasheet

Rtf1 Antibody (orb1530582)

Description	Rtf1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 660 and C-terminus
Target	Rtf1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q92541
NCBI	NP_055953.3
Dilution Range	WB - 1:10,000 - 1:30,000; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human Rtf1 by western blot....



Detection of human Rtf1 by western blot ...