



Product Datasheet

SMC4 Antibody (orb1530738)

Description	SMC4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Between 600 and 650
Target	SMC4
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9NTJ3
NCBI	NP_005487.2
Dilution Range	WB - 1:1,000 - 1:10,000; ICC-IF - 1:500 - 1:2,000
Expiration Date	12 months from date of receipt.



Detection of human SMC4 by western blot....



Localization of SMC4 to Mitotic Chromoso...