


Product Datasheet

SMC3, Phospho (S383) Antibody (orb1530740)

Description	SMC3, Phospho (S383) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Surrounding serine 383
Target	SMC3, Phospho (S383)
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9UQE7
NCBI	NP_005436.1
Dilution Range	WB - Not recommended; IP - 2 - 6 µg/mg lysate
Expiration Date	12 months from date of receipt.



Immunoprecipitation
of human Phospho
SMC3...