


Product Datasheet

XRCC1, Phospho (S518/T519/T523) Antibody (orb1530744)

Description	XRCC1, Phospho (S518/T519/T523) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Surrounding serine 518 and threonine 519 /523
Target	XRCC1, Phospho (S518/T519/T523)
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P18887
NCBI	NP_006288.2
Dilution Range	WB - 1:1,000 - 1:5,000; IHC - 1:2,000 - 1:10,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human Phospho-ERCC1 (S518/T519/T523)