


Product Datasheet

BRCA1, Phospho (S1524) Antibody (orb1530810)

Description	BRCA1, Phospho (S1524) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Surrounding serine 1524
Target	BRCA1, Phospho (S1524)
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P38398
NCBI	NP_009225.1
Dilution Range	WB - 1:500 - 1:5,000
Expiration Date	12 months from date of receipt.



Detection of
human
Phospho
BRCA1
(S1524)...