

## Product Datasheet

### E Tag Antibody (orb1530845)

<b>Description</b>	E Tag Antibody
<b>Species/Host</b>	Goat
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	ICC, IF
<b>Target</b>	E Tag
<b>Preservatives</b>	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
<b>Form/Appearance</b>	Whole IgG
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	2 - 8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Format: Whole IgG
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	IHC - Not applicable; ICC - 1:100 - 1:400; IF - 1:100 - 1:400
<b>Expiration Date</b>	12 months from date of receipt.