


Product Datasheet

T7 Tag Antibody (orb1530864)

Description	T7 Tag Antibody
Species/Host	Goat
Conjugation	HRP
Tested Applications	ELISA, ICC, WB
Target	T7 Tag
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.05% Pro-Clean 400
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	WB - 1:1,000 - 1:30,000; IHC - Not applicable; ICC - 1:200 - 1:500; ELISA - 1:10,000 - 1:100,000
Expiration Date	12 months from date of receipt.



Detection of T7-tagged protein by wester...