


Product Datasheet

6-His Tag Antibody (orb1530868)

Description	6-His Tag Antibody
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	ELISA, WB
Target	6-His Tag
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	WB - 1:10,000 - 1:100,000; ELISA - 1:10,000 - 1:100,000
Expiration Date	12 months from date of receipt.



Detection of 6-His-tagged protein by WB....