

Product Datasheet

6-His Tag Antibody (orb1530871)

Description	6-His Tag Antibody
Species/Host	Goat
Conjugation	FITC
Tested Applications	ICC, IF
Target	6-His Tag
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	IHC - Not applicable; ICC - 1:100 - 1:400; IF - 1:100 - 1:400
Expiration Date	12 months from date of receipt.