


Product Datasheet

Myc Tag Antibody (orb1530877)

Description	Myc Tag Antibody
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	ELISA, ICC, WB
Target	Myc Tag
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.05% Pro-Clean 400
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P01106
Dilution Range	WB - 1:1,000 - 1:30,000; IHC - Not applicable; ICC - 1:200 - 1:500; ELISA - 1:10,000 - 1:100,000
Expiration Date	12 months from date of receipt.



Detection of c-myc-tagged protein by WB....